

5/28/2021

***C.elegans* mtDNA measurement by qPCR**

by Qian Zhang

Ye Tian Lab, IGDB, CAS.

Preparation of DNA template

1. 30 early L4 worms were collected to a PCR tube containing 50 μ l WLB.
2. Incubate 100 min at 65°C.
3. Heat inactivate PK at 95°C for 15min.
4. Cool on ice.

Setup of qPCR reaction

5. Setup reaction for each column:
 - primer mix, 4 μ l
 - SYBR Green Real-time PCR Master Mix (TOYOBO QPK-201), 10 μ l
 - worm lysate, 3-4 μ l
 - Add ddH₂O to final volume of 20 μ l
6. Both mtDNA and nDNA needed to be measured in each lysate for 3-4 technical repeats.
7. Primers:

mtDNA-F	GCTTTTCTTTATATGTTTTGTG
mtDNA-R	TCACCTTCAGAAAAATCAAATGG
nDNA-F	AGGCTAAGCCGGGGTAAGTT
nDNA-R	GCCAAAAGCTTAAACTGCGG
8. The mtDNA/nDNA ratio was then calculated using the comparative $\Delta\Delta$ Ct method.

Worm Lysis Buffer (WLB):

50 mM KCl,
10 mM Tris-HCl (pH 8.3),
2.5 mM MgCl₂,
0.45% NP-40,
0.45% Tween 20,
0.01% gelatin,
fresh proteinase K (200 μ g/ml)